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BTX[®]

a division of Harvard Bioscience, Inc.

NOTE: Products in this catalog are for Research Use Only. Not for use on humans unless proper investigational device regulations have been followed.

Harvard is a registered trademark of Harvard University. The marks Harvard Apparatus and Harvard Bioscience are being used pursuant to a license agreement between Harvard University and Harvard Bioscience, Inc.

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
Typical Results




Figure 1: Field strength and time constants for gram-negative bacteria in 1 mm gap cuvettes using the ECM 630 and 399 models. Field strengths and time constants remain grouped around discrete values. The ECM 630 would be a good choice for labs that are currently doing simple transformations now, but plan on working with other cell lines in the future, while the ECM 399 would be ideal for the lab that is only interested in simple transformations of bacteria and yeast where the pulse duration is no longer than 5 - 6 ms.

Figures 2 and 3: Field Strength and pulse duration values for gram-positive bacteria in 1 mm and 2 mm cuvettes. Unlike the gram-negative bacteria the field strengths and time constants are more variable with gram-positive. The ECM 630 is flexible and settings can be adjusted for optimizing multiple cell lines.

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References

- Cashmore, TJ, et al. Identification of a Membrane Protein Required for Lipomannan Maturation and Lipoarabinomannan Synthesis in Corynebacterineae. *Biol. Chem.* 2017; jbc-M116.
- Lee JH, Park SW, Kim YM, Oh JI. Functional Characterization of the CutI Gene for the Transcription of Carbon Monoxide Dehydrogenase Genes in Mycobacterium Sp. Strain JC1 DSM 3803. *J Microbiol.* 2017;55: 31-36.
- Weger-Lucarelli, J, et al. Development and Characterization of Recombinant Virus Generated from a New World Zika Virus Infectious Clone. *J. Virol.* 2017 Oct;91(19): e01765-16.
- Bakhshi PK, Bain J, Gul MO, Stride E, Edirisinghe M, Staniland SS. Manufacturing Man - Made Magnetosomes: High - Throughput In Situ Synthesis of Biomimetic Magnetite Loaded Nanovesicles. *Macromol Biosci.* 2016;16: 1555-1561.
- Kim, SJ, Ohgwe K, Cerniglia, CE. Phenotype-Based Identification of Key Enzymes for Polycyclic Aromatic Hydrocarbon (PAH) Metabolism from Mycobacteria Using Transposon Mutagenesis and a PAH Spray Plate. In: McGenity T., Timmis K., Nogales B. (eds) *Hydrocarbon and Lipid Microbiology Protocols.* Springer Protocols Handbooks. 2015. Springer, Berlin, Heidelberg. 161-175.
- Lee J, Jang YS, Papoutsakis ET, Lee SY. Stable and Enhanced Gene Expression in *Clostridium Acetobutylicum* using Synthetic Untranslated Regions with a Stem-Loop. *J Biotech.* 2016; 230: 40-43.
- Lin TL, Pan YJ, Hsieh PF, Hsu CR, Wu MC, Wang JT. Impenem Represses CRISPR-Cas Interference of DNA Acquisition through H-NS Stimulation in *Klebsiella Pneumoniae*. *Sci Rep.* 2016; 17: 21644.
- Oakes BL, et al. Profiling of Engineering Hotspots Identifies an Allosteric CRISPR-Cas9 Switch. *Nat Biotechnol* 2016; 34: 646-651.
- Rojas-Sánchez S, et al. Transcription of *Leishmania major* U2 small nuclear RNA gene is directed by extragenic sequences located within a tRNA-like and a tRNA-Ala gene. *Parasite Vector* 2016; 9: 401.
- Wisniewski JA, et al. TcpM: a novel relaxase that mediates transfer of large conjugative plasmids from *Clostridium perfringens*. *Mol. Microbiol* 2016 Mar;99(5): 884-96.
- Dobroff AS, et al. Ligand - Directed Profiling of Organelles with Internalizing Phage Libraries. *Curr Protoc Protein Sci.* 2015;79: 30.4.1-30.4.30.
- Jones SC, Price CT, Santic M, Abu Kwaik Y. Selective Requirement of the Shikimate Pathway of *Legionella Pneumophila* for Intravacuolar Growth within Human Macrophages but not with in *Acanthamoeba*. *Infect Immun* 2015;83: 2487-2495.
- Kim S, et al. Redox-Switch Regulatory Mechanism of Thiolase from *Clostridium Acetobutylicum*. *Nat Commun* 2015;22 8410.
- Knut M, Tolstorkov I, Cregg J. Electroporation of *Pichia pastoris*. In: van den Berg M., Maruthachalam K. (eds) *Genetic Transformation Systems in Fungi Volume 1.* Fungal Biology. 2015. Springer, Cham 87-91.
- Ren J, et al. Acetylation Regulates Survival of *Salmonella Enterica* Serovar Typhimurium under Acid Stress. *Appl Environ Microb* 2015;81(17): 5675-5682.
- Wisniewski JA, Teng WL, Bannam TL, Rood JI. Two Novel Membrane Proteins, TcpD and TcpE, Are Essential for Conjugative Transfer of pCW3 in *Clostridium Perfringens*. *J Bacteriol.* 2015;197: 774-781.
- Cameron DE, Collins JJ. Tunable Protein Degradation in Bacteria. *Nat Biotechnol* 2014;32: 1276-1281.
- Ravi Kant H, Balamurali M, Meenakshisundaram S. Enhancing precursors availability in *Pichia pastoris* for the overproduction of S-adenosyl-methionine employing molecular strategies with process tuning. *J. Biotech.* 2014;188: 112-121.
- Nováková J, et al. Improved method for high-efficiency electrotransformation of *Escherichia coli* with the large BAC plasmids. *Folia Microbiol* 2014; 59(1): 53-61.
- Ueno Y, Hoshino K, Tagawa Y. Mutations in the major outer membrane protein gene from *Histophilus somni* by an allelic exchange method. *J. Microbiol. Meth.* 2014;106: 83-92.
- Wang L, et al. A new strategy for secretory expression and mixed fermentation of recombinant human collagen $\alpha 1$ (III) chain in *Pichia pastoris*. *Biotech. Bioproc. Eng.* 2014;19(5): 916-924.
- Wei Y, Hu H, Lun ZR, Li Z. Centrin3 in Trypanosomes Maintains the Stability of a Flagellar Inner-Arm Dynein for Cell Motility. *Nat Comms.* 2014; 5: 4060.

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
Typical Results

Electroporation of anthers (top row) prior to their culture induces faster growth (c), somatic embryo formation (e) and, ultimately, haploid plant regeneration (g), as shown here for field pea.

Ravinder Kaur Grewal, Monika Lulsdorf, Janine Croser, Sergio Ochatt, Albert Vanderberg, Thoma D. Warkentin., Doubled-haploid Production in Chickpea *Cicer arietinum* L.): Role of Stress Treatments. *Plant Cell Reports*, 30, May, 2009

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References

Ruiz N, et al. Non-Invasive Delivery of dsRNA into De-Waxed Tick Eggs by Electroporation. PLoS ONE . 2015; 10:e0130008.

Wojcik A, Rybczynski JJ. Electroporation and morphogenic potential of Gentiana kurroo (Royle) embryogenic cell suspension protoplasts. Journal of Biotech Computat Biol Bionanotech 2015;96(1): 19-29.

Gorchakov RV, Tesh RB, Weaver SC, Nasar F. Generation of an infectious Negev virus cDNA clone. J Gen Virol. 2014; 95(9):2071-2074.

Jeon K, Suresh A, Kim Y C. Highly efficient molecular delivery into Chlamydomonas reinhardtii by electroporation Korean J Chem Eng . 2013;30(8): 1626-1630.

Suresh A, Kim YC. Translocation of cell penetrating peptides on Chlamydomonas reinhardtii. Biotechnol Bioeng 2013;110(10): 2795-2801.

Karim S, Troiano E, Mather TN. Functional genomics tool: gene silencing in Ixodes scapularis eggs and nymphs by electroporated dsRNA. BMC Biotechnol, 2010;10(1): 1.

Thomas, J. Electroporation, an Alternative to Biolistics for Transfection of Bombyx mori Embryos and Larval Tissues. J Insect Sci . 2003; 3:17.

Ravinder Kaur Grewal, Monika Lulsdorf, Janine Croser, Sergio Ochatt, Albert Vanderberg, Thoma D. Warkentin., Doubled-haploid Production in Chickpea Cicer arietinum L.): Role of Stress Treatments. Plant Cell Reports, 30, May, 2009

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

Typical Results

PC12 cells were transfected with eGFP plasmid DNA using the ECM 830 system. GFP expression was observed 48 hrs. post transfection. Courtesy of Anne Chiamello, Ph.D Associate Professor, George Washington University Medical Center.

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References

- Bongard N, et al. Interference of Retroviral Envelope with Vaccine-Induced CD8+ T Cell Responses is Relieved by Co-Administration of Cytokine-Encoding Vectors. *Retrovirology* 2017;14: 28.
- Callif BL, Maunze B, Krueger NL, Simpson MT, Blackmore MG. The application of CRISPR technology to high content screening in primary neurons. *Mol Cell Neurosci*. 2017;80: 170-179.
- Denies S, et al. In Vitro Exploration of a Myeloid-Derived Suppressor Cell Line as Vehicle for Cancer Gene Therapy. *Cancer Gene Ther* . 2017;24: 149-155.
- Eyquem J, et al. Targeting a CAR to the TRAC Locus with CRISPR/Cas9 Enhances Tumour Rejection. *Nature*. 2017;543: 113-117.
- Huang J, et al. Generation and Comparison of CRISPR-Cas9 and Cre-Mediated Genetically Engineered Mouse Models of Sarcoma. *Nature Comm*. 2017;8: 15999.
- Kuo CY, Lin CH, Hsu T. VHL inactivation in precancerous kidney cells induces an inflammatory response via ER stress-activated IRE1 α signaling. *Cancer Res* . 2017;77: 3406-3416.
- Liu, X et al. Novel T cells with improved in vivo anti-tumor activity generated by RNA electroporation. *Protein Cell* 2017 Jul;8(7): 514-526.
- Lopes A, Vanvarenberg K, Preat V, Vandermeulen G. Codon-Optimized P1A-Encoding DNA Vaccine: Toward a Therapeutic Vaccination against P815 Mastocytoma. *Mol Ther Nucleic Acids*. 2017 Sep 15; 8: 404-415.
- Ren J, et al. A versatile system for rapid multiplex genome-edited CAR T cell generation. *Oncotarget*. 2017;8(10): 17002.
- Seyid Jafari SM, et al. Efficacy of In Vivo Electroporation-Mediated IL-10 Gene Delivery on Survival of Skin Flaps. *J Membrane Biol*. 2017; 10.1007/s00232-017-9974-x.
- Sun Y, et al. Coinjection of IL2 DNA Enhances E7-Specific Antitumor Immunity Elicited by Intravaginal Therapeutic HPV DNA Vaccination with Electroporation. *Gene Ther*. 2017;24: 408-415.
- Xie Z, et al. Optimization of a CRISPR/Cas9-mediated Knock-in Strategy at the Porcine Rosa26 Locus in Porcine Foetal Fibroblasts. *Sci Rep* . 2017;7: 02785-y.

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- An L, et al. Efficient generation of FVII gene knockout mice using CRISPR/Cas9 nuclease and truncated guided RNAs. *Sci Rep* . 2016;6: 25199.
- Davtyan H, et al. BTX AgilePulse System is an Effective Electroporation Device for Intramuscular and Intradermal Delivery of DNA Vaccine. *Sci Rep* . 2016;6: 28912.
- Hawksworth DJ, Ziemann R, Manoj S. Enhanced Immune Response Following DNA Immunization Through In Vivo Electroporation. *Monoclon Antib Immunodiagn Immunother* 2016;35(5): 239-244.
- Hoppe PS, et al. Early Myeloid Lineage Choice is not Initiated by Random PU.1 to GATA1 Protein Ratios. *Nature*. 2016 Jul 14;535(7611): 299-302.
- Martell JD, et al. A Split Horseradish Peroxidase for the Detection of Intercellular Protein – Protein Interactions and Sensitive Visualization of Synapses. *Nature Biotech*. 2016 Jul;34(7): 774-780.
- Peterson CW, et al. Long-Term Multi-Lineage Engraftment of Genome-Edited Hematopoietic Stem Cells after Autologous Transplantation in Nonhuman Primates. *Blood*. 2016; 127: 2416-2426.
- Stanga S, et al. APP-dependent glial cell line-derived neurotrophic factor gene expression drives neuromuscular junction formation. *FASEB J* . 2016 May;30(5), 1696-1711.
- Wang K, et al. Efficient Generation of Orthologous Point Mutations in Pigs via CRISPR-assisted ssODN-mediated Homology-Directed Repair. *Mol Ther Nucleic Acids* . 2016 Nov;5(11): e396.
- Yamasaki S, et al. In Vivo Dendritic Cell Targeting Cellular Vaccine Induces CD4+ Tfh Cell-Dependent Antibody against Influenza Virus. *Sci Rep* . 2016; 6: 35173.
- Kim NY, et al. Enhanced Immune Response to DNA Vaccine Encoding Bacillus Anthracis PA-D4 Protects Mice Against Anthrax Spore Challenge. *PLoS ONE* . 2015 Oct 2;10(10): e0139671.
- Wang J, et al. Homology-Driven Genome Editing in Hematopoietic Stem and Progenitor Cells using ZFN mRNA and AAV6 Donors. *Nat Biotechnol* 2015 Dec;33(12): 1256-1263.

Please note: Nucleofector is a registered trademark of Lonza.

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Typical Results


Figure 1: Bioluminescence Imaging of Tibialis anterior muscles injected with a luciferase plasmid and electroporated compared to injection alone. Provided by Carmen Bertoni Ph.D, Department of Neurobiology UCLA.

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
Figure 2: Stable long term expression obtained after intramuscular injection of a PhiC31 plasmid encoding the green fluorescent protein (GFP) under the control of a muscle specific promoter followed by electroporation. Provided by Carmen Bertoni Ph.D Dept. of Neurobiology UCLA.

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
Figure 3: A. Sagittal views of a P15 mouse 2 weeks post electroporation; 60 μ m sections. In addition to the RMS, GFP-positive cells can be found throughout the olfactory bulb. B/C. Enlarged views of the boxed regions in (a) showing olfactory bulb (b), and the RMS (c). Inset in (b) highlights a periglomerular cell. (Chesler et al. 2008).

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References

- Chu SF, et al. Upregulating the expression of survivin-HBXIP complex contributes to the protective role of IMM-H004 in transient global cerebral ischemia/reperfusion. *Mol Neurobiol.* 2017;54(1): 524-540.
- Edelblute C M, et al. Controllable Moderate Heating Enhances the Therapeutic Efficacy of Irreversible Electroporation for Pancreatic Cancer. *Sci Rep.* 2017;7: 11767.
- Fekete CD, et al. In vivo transgenic expression of collybistin in neurons of the rat cerebral cortex. *J Comp Neurol.* 2017 Apr 1;525(5): 1291-1311.
- Homman-Ludiyé J, Kwan WC, de Souza MJ, Rodger J, Bourne J.A. Ephrin-A2 regulates excitatory neuron differentiation and interneuron migration in the developing neocortex. *Neurosci.* 2017;7: 11813.
- Ji L, et al. Defective neuronal migration and inhibition of bipolar to multipolar transition of migrating neural cells by Mesoderm-Specific Transcript, Mest, in the developing mouse neocortex. *Neurosci.* 2017;355: 126-140.
- Zhou Z, et al. Rab28 is a TBC1D1/TBC1D4 substrate involved in GLUT4 trafficking. *FEBS Lett.* 2017;591(1): 88-96.
- Gordy JT, Luo K, Zhang H, Biragyn A, Markham RB. Fusion of the dendritic cell-targeting chemokine MIP3 α to melanoma antigen Gp100 in a therapeutic DNA vaccine significantly enhances immunogenicity and survival in a mouse melanoma model. *J Immunother Cancer.* 2016;4(1): 96.
- Ji L, Kim NH, Huh SO, Rhee HJ. Depletion of Inositol Polyphosphate 4-Phosphatase II Suppresses Callosal Axon Formation in the Developing Mice. *Mol Cells.* 2016;39(6): 501-507.
- Pacary E, Guillemot F. Cerebral cortex electroporation to study projection neuron migration. *Curr Protoc Neurosci.* 2016 Oct 3;77: 1-2.26. 18.
- Stanga S, et al. APP-dependent glial cell line-derived neurotrophic factor gene expression drives neuromuscular junction formation. *EASEB J.* 2016;30(5): 1696-1711.
- Varela JA. Targeting neurotransmitter receptors with nanoparticles in vivo allows single-molecule tracking in acute brain slices. *Nat Comm.* 2016;7: 10947.
- Usmani A, et al. A non-surgical approach for male germ cell mediated gene transmission through transgenesis. *Sci Rep.* 2013;3: 3430.
- Dmitriev RI, et al. Imaging oxygen in neural cell and tissue models by means of anionic cell-permeable phosphorescent nanoparticles. *Cell Mol Life Sci.* 2015;72(2): 367-381.
- Peng H, et al. Nlrp2, a Maternal Effect Gene Required for Early Embryonic Development in the Mouse. *PLoS One.* 2012;7: e30344.
- Xu H, et al. The Calmodulin-Stimulated Adenylate Cyclase ADCY8 Sets the Sensitivity of Zebrafish Retinal Axons to Midline Repellents and Is Required for Normal Midline Crossing. *J Neurosci.* 2010;30: 7423-7433.
- Chesler, et al. Selective Gene Expression by Postnatal Electroporation during Olfactory Interneuron Neurogenesis. *PLoS One.* 2008;3(1): e1517
- Tonelli D, et al. Single-Cell Detection of microRNAs in Developing Vertebrate Embryos After Acute Administration of a Dual-Fluorescence Reporter/Sensor Plasmid. *BioTechniques.* 2006 Dec;41: 727-732.

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

Typical Results

Fusion of dendritic cells to tumor cells was determined by loading either CMFDA or CMTMR fluorescent labels into cells. Fusion was carried out by using PEG-mediated fusion and Electrofusion methods. The cells were analysed by flow cytometry. Cells expressing both fluorescent labels were considered fused. Electrofusion resulted in significantly higher numbers of fused cells in comparison to PEG.

Rimas J. Orentas, Dennis Schauer, Qian Bin, and Bryon D. Johnson; Electrofusion of a weakly immunogenic Neuroblastoma with Dendritic Cells Produces a Tumor Vaccine. *Cell. Immunol* 2001;213: 4-13.

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Sper RB, et al. Generation of a Stable Transgenic Swine Model Expressing a Porcine Histone 2B-eGFP Fusion Protein for Cell Tracking and Chromosome Dynamics Studies. *PLoS ONE* . 2017; 12(1): e0169242.

Huan YJ, et al. Alteration of the DNA methylation status of donor cells impairs the developmental competence of porcine cloned embryos. *J Reproduc Develop* . 2016; 62(1): 71-77.

Jeong YH, et al. Knock-in fibroblasts and transgenic blastocysts for expression of human FGF2 in the bovine β -casein gene locus using CRISPR/Cas9 nuclease-mediated homologous recombination. *Zygote* . 2016; 24(3): 442-456.

Sun Y, et al. Deletion of a Yci1 domain protein of *Candida albicans* allows homothallic mating in MTL heterozygous cells. *mBio*. 2016;7(2): e00465-16.

Thornburg NJ, et al. H7N9 influenza virus neutralizing antibodies that possess few somatic mutations. *J Clin Invest*. 2016; 126(4): 1482-1494.

Wang K, et al. Efficient generation of orthologous point mutations in pigs via CRISPR-assisted ssODN-mediated homology-directed repair. *Mol Ther Nucleic Acids* . 2016 Nov;5(11): e396.

Mangan PR, et al. Dual inhibition of interleukin-23 and interleukin-17 offers superior efficacy in mouse models of autoimmunity. *J Pharmacol Exp Ther* , 2015; 354(2): 152-165.

Wang K, et al. Efficient generation of myostatin mutations in pigs using the CRISPR/Cas9 system. *Sci Reports* . 2015; 5, 16623.


Zhang Y, Luo W, Wang Y, Liu Y, Zheng L. Purified dendritic cell-tumor fusion hybrids supplemented with non-adherent dendritic cells fraction are superior activators of antitumor immunity. *PLoS ONE* . 2014;9(1): e86772.

Tan C, et al (2013). Impact of anti-CD25 monoclonal antibody on dendritic cell-tumor fusion vaccine efficacy in a murine melanoma model. *J Transl Med* . 2013; 11(1): 148.

Parkhurst MR, et al. Hybrids of Dendritic Cells and Tumor Cells Generated by Electrofusion Simultaneously Present Immunodominant Epitopes from Multiple Human Tumor-Associated Antigens in the Context of MHC Class I and Class II Molecules. *J Immunol* 2003;170: 5317-5325.


Orentas R, et al. Electrofusion of a Weakly Immunogenic Neuroblastoma with Dendritic Cells Produces a Tumor Vaccine. *Cell Immunol* 2001;213: 4-13.

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Typical Results

Simple direct scale-up of transfection using the AgilePulse MAX System. K562 cells were transfected with GFP+ mRNA in standard 4 mm gap electroporation laboratory cuvettes (0.5 ml) and AgilePulse MAX chambers (6.0 ml). Identical pulse parameters were applied to all three electroporation chambers. At 24 hours post electroporation, the percent transfection was determined by flow cytometry. The efficiency of transfection was comparable for both chambers.

(Markovic S, et al. Preparing clinical-grade myeloid dendritic cells by electroporation-mediated transfection of in vitro amplified tumor-derived mRNA and safety testing in stage IV malignant melanoma. J Translat Med. 2006 Aug;4: 35.)

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
Typical Results

Dermal delivery of DNA vaccine with electroporation using the AgilePulse In Vivo System increases the immune response of PSA specific CD8+ T cells over intradermal (ID) alone or intramuscular injection (IM), with or without electroporation, (Roos, A, et al. Mol. Ther. 2006;13(2): 320-327.)

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
Typical Results

Use of ECM 830 electroporation system and Tweezertrode electrodes for in vivo genome editing via microhomology-mediated end-joining (MMEJ)-mediated targeted integration. (a) Experimental scheme for targeted Actb-2A-mCherry knock-in in fetal brain via in utero electroporation. (b) Representative immunofluorescence images of neurons showing correct mCherry knock-in at the Actb locus. Scale bar, 50 μ m. GFP, transfected cells. (c) Relative knock-in efficiency measured by the percentage of mCherry⁺ cells among GFP⁺ cells. Adapted from EBioMedicine, Xuan Yao, et al. (2017) with permission from Elsevier.

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References

- Yao X, et al. CRISPR/Cas9 - Mediated Precise Targeted Integration In Vivo Using a Double Cut Donor with Short Homology Arms. EBioMedicine. 2017; 20:19-26. This reference must be listed here or as part of the figure caption.
- Qin W, et al. Efficient CRISPR/Cas9-Mediated Genome Editing in Mice by Zygote Electroporation of Nuclease. Genetics. 2015; 200:423-430.
- Sidik S, et al. Efficient Genome Engineering of Toxoplasma gondii Using CRISPR/Cas9. PLoS ONE. 2014; 9: e100450.
- Wang X, et al. Efficient CRISPR/Cas9-Mediated Biallelic Gene Disruption and Site-Specific Knockin after Rapid Selection of Highly Active sgRNAs in Pigs. Sci Rep. 2015; 5:13348.
- Xie Z, et al. Optimization of a CRISPR/Cas9-Mediated Knock-In Strategy at the Porcine Rosa26 Locus in Porcine Foetal Fibroblasts. Sci Rep. 2017; 7:3036.
- Dimitrov L, et al. Germline Gene Editing in Chickens by Efficient CRISPR-Mediated Homologous Recombination in Primordial Germ Cells. PLoS ONE. 2016; 11: e0154303.
- Shin SE, et al. CRISPR/Cas9-Induced Knockout and Knock-in Mutations in Chlamydomonas reinhardtii. Sci Rep. 2016; 6:27810.
- An L, et al. Efficient Generation of FVII Gene Knockout Mice using CRISPR/Cas9 Nuclease and Truncated Guided RNAs. Sci Rep. 2016; 6:25199.

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Four different transgenic mice expressing Abs to human Ag were used to compare the efficiency of electrofusion to PEG fusion. Transgenic, human-Ab producing mice immunized with tetanus toxoid (TT) provided spleens for fusion to SP2/O mouse myeloma cells. PEG fusion was performed with standard protocols. For electrofusion, mouse and SP2/O cells were washed twice in Cytofusion medium then mixed in the fusion chamber and tri-phasic pulse applied. Cells were recovered after 30 min and cultured in 96-well plates at 5000 cells/ml. Antigen-specific clones were counted using ELISA or HTRF, normalized to 100M cells. Source: Data courtesy of M. Coccia, PhD, Platform Development Group, Medarex Inc, Milpitas, CA)

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
Electroporation increases the immune response. C57Bl/6 mice were immunized once with 10 pVax-PSA/20 ml PBS intradermally (id) on each flank with or without electroporation (EP) or intramuscular (IM) in each TA muscle. Blood was collected on days 11, 13 and 15 and effector cells were stimulated for 4 h with 100 nM PSA-derived peptide psa65-73 or a control peptide GP33. The activated CD8+ T cells were quantified by intracellular cytokine staining for IFN- γ and analyzed by flow cytometry. The data demonstrate that dermal delivery of DNA vaccine along with electroporation using the AgilePulse In Vivo System increases the immune response of PSA-specific CD8+ T cells over intradermal (ID) alone or intramuscular injection (IM), with or without electroporation. Roos, A, et al. *Molec. Ther.* 2006;13(2): 320-327.

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货号	描述	包含内容
47-0400N	AgilePulse In Vivo ID (Intra-Dermal) 系统	AgilePulse ID电穿孔仪, 平行针形矩阵电极手柄 (47-0000), 4针阵列电极 6针阵列电极和操作软件
47-0401N	AgilePulse In Vivo ID (Intra-Dermal)电穿孔仪	仅含AgilePulse In Vivo ID (Intra-Dermal) 电穿孔仪
47-0500N	AgilePulse In Vivo IM (Intra-Muscular) 系统	AgilePulse IM 电穿孔仪, 平行针形阵列电极手柄 (47-0000), 4-针阵列电极 6针阵列电极和操作软件
47-0501N	AgilePulse In Vivo IM (Intra-Muscular)电穿孔仪	仅含AgilePulse In Vivo IM (Intra-Muscular) 电穿孔仪
47-0090	AgilePulse In Vivo电极适配器	AAgilePulse In Vivo电极适配器

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
Efficiency of GFP RNA Delivery is the same in 4 mm Cuvettes (0.5 ml) and AP Max Chamber (5 ml)

Scale-up of K562 cell mRNA transfection using AgilePulse MAX system. K562 cells (myelomonocytic cells commonly used for natural killer target cell assays) were transfected with GFP+ mRNA in both laboratory cuvettes (0.5 mL) and the large volume AgilePulse MAX System. Cells were suspended to a cell density of 20 million cells/mL in Cytoporation Medium. GFP mRNA was added to the cell suspension to a final concentration of 40 µg/mL. The small 0.5 mL volume transfections were carried out using standard 4 mm gap electroporation cuvettes. The large volume transfections were performed with the AgilePulse MAX using a 4 mm gap large volume chamber. Identical pulse parameters were applied to all three volumes. At 24 hours post-electroporation, the percent transfection was determined by flow cytometry. The efficiency of transfection was comparable for all volumes. (Markovic S, et al. Preparing clinical-grade myeloid dendritic cells by electroporation-mediated transfection of in vitro amplified tumor-derived mRNA and safety testing in stage IV malignant melanoma. J Translat Med. 2006 Aug;4: 35.)

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