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Alamar Blue Cell Viability Reagent

Introduction

Cell health can be monitored by numerous methods. Plasma membrane integrity, DNA synthesis, DNA content, enzyme activity, presence of ATP, and cellular reducing conditions are known indicators of cell viability and cell death. Alamar Blue cell viability reagent functions as a cell health indicator by using the reducing power of living cells to quantitatively measure the proliferation of various human and animal cell lines, bacteria, plant, and fungi allowing you to establish relative cytotoxicity of agents within various chemical classes. When cells are alive they maintain a reducing environment within the cytosol of the cell. Resazurin, the active ingredient of Alamar Blue reagent, is a non-toxic, cell permeable compound that is blue in color and virtually non-fluorescent. Upon entering cells, resazurin is reduced to resorufin, a compound that is red in color and highly fluorescent. Viable cells continuously convert resazurin to resorufin, increasing the overall fluorescence and color of the media surrounding cells.

Alamar Blue cell viability reagent is used to assess cell viability by simply adding the 10X, ready-to-use solution to mammalian or bacterial cells in culture media. There is no requirement to aspirate media from cells or place cells in minimal media. Consequently, Alamar Blue reagent can easily be used in a single tube or microtiter plate format in a “no wash” fashion. Simply add Alamar Blue reagent as 10% of the sample volume (i.e., add 10 μ L Alamar Blue reagent to 100 μ L sample), followed by a 1–4 hours incubation at 37°C. Longer incubation times may be used for greater sensitivity without compromising cell health. The resulting fluorescence is read on a plate reader or fluorescence spectrophotometer. Alternatively, the absorbance of Alamar Blue reagent can be read on a spectrophotometer. Finally, results are analyzed by plotting fluorescence intensity (or absorbance) versus compound concentration.

Materials required but not provided

Mammalian or bacterial cells in appropriate medium
Appropriate 96- or 384-well plates
Optional: 3% SDS in phosphate buffered saline (PBS), pH 7.4

Preparing Cells

Mammalian Cells—Adherent: Plate mammalian cells in a cell culture flask or dish, and allow cells to adhere and grow for approximately 4–24 hours at 37°C and 5% CO₂ before proceeding with the assay.

Mammalian Cells—Suspension: Plate mammalian cells in a cell culture flask or dish, and use cells immediately for the assay or allow cells to grow for up to 24 hours at 37°C and 5% CO₂ before proceeding with the assay.

Bacterial Cells: For details, see references 2 and 3. Notes Alamar Blue reagent is stable to multiple freeze/thaw cycles and its activity is not affected if the reagent is frozen.

General Guidelines

Cell types assayed with Alamar Blue reagent include mammalian, bacterial (including biofilms), plant, and fish cells. More specifically Alamar Blue reagent has been tested on hepatocytes, such as HepG2 cells, as well as cells of primary origin.

Be sure to include appropriate assay controls. To minimize experimental errors, we recommend making measurements from a minimum of 4–8 replicates of experimental and no-cell control samples.

You may need to determine the plating density and incubation time for the Alamar Blue assay for each cell type and use conditions such that the assay is in the linear range.

If you plan to use longer incubation time (overnight), be sure to maintain sterile conditions during reagent addition and incubation to avoid microbial contaminants. Contaminated cultures will yield erroneous results as microbial contaminants also reduce Alamar Blue reagent.

Fetal bovine serum (FBS) and bovine serum albumin (BSA) cause some quenching of fluorescence. We recommend using the same serum concentration in controls to account for this quenching. Other media components, such as phenol red do not interfere with the assay.

Protocol

Optional: Treat cells with the test compound 24–72 hours prior to performing the Alamar Blue cytotoxicity assay. Add 1/10th volume of Alamar Blue reagent directly to cells in culture medium.

Incubate for 1 to 4 hours at 37°C in a cell culture incubator, **protected from direct light**.

Note: Sensitivity of detection increases with longer incubation times. For samples with fewer cells, use longer incubation times of up to 24 hours.

Record results using fluorescence or absorbance as follows:

Fluorescence: Read fluorescence using a fluorescence excitation wavelength of 540–570 nm (peak excitation is 570 nm). Read fluorescence emission at 580–610 nm (peak emission is 585 nm).

Absorbance: Monitor the absorbance of Alamar Blue at 570 nm, using 600 nm as a reference wavelength (normalized to the 600 nm value).

Note: Fluorescence mode measurements are more sensitive. When fluorescence instrumentation is unavailable, monitor the absorbance of Alamar Blue reagent. Assay plates or tubes can be wrapped in foil, stored at 4°C, and read within 1–3 days without affecting the fluorescence or absorbance values.

Optional: Add 50 μ L 3% SDS directly to 100 μ L of cells in Alamar Blue reagent to stop the reaction.

References

1. Invest Ophthalmol Vis Sci 38, 1929 (1997); 2. Infect Immun 65, 3193 (1997); 3. J Antimicrob Chemother 57, 1100 (2006); 4. Phytochem Anal 12, 340 (2001); 5. Anal Biochem 344, 76 (2005)
FAQ's

General Questions

Q: How does Alamar Blue work?

A: Healthy living cells maintain a reducing state within their cytosol. This “reducing potential” of cells converts Alamar Blue reagent into a detectable fluorescent (or absorbent) product.

Q: Is Alamar Blue reagent toxic?

A: No. Alamar Blue reagent is a safe, non-toxic reagent to both the sample and user.

Q: Does Alamar Blue reagent need reconstitution?

A: No, Alamar Blue reagent is supplied as a 10X, ready-to-use solution.

Q: Can I use Alamar Blue reagent with suspension cells too?

A: Yes. Alamar Blue reagent works on adherent and suspension mammalian cells.

Q: Can I use Alamar Blue reagent with non-mammalian cells, such as bacteria?

A: Yes, Alamar Blue reagent has been shown to work with bacterial and plant cells.

Q: Alamar Blue reagent is not the most expensive cytotoxicity indicator on the market, does that mean it doesn't work as well as other reagents?

A: Actually, Alamar Blue reagent is comparable to other often more expensive cytotoxicity indicators.

Q: Since Alamar Blue is an absorbance or fluorescence readout, is it as sensitive as a luminescence product?

A: Alamar Blue reagent is sensitive enough to detect less than 50 mammalian cells in a single well of a 96-well plate.

Storage Questions

Q: What if I left the Alamar Blue stock reagent at room temperature, overnight?

A: The reagent is stable for up to 12 months when stored at room temperature (~22°C).

Q: I accidentally froze the Alamar Blue stock reagent, can I still use it?

A: Yes. Alamar Blue reagent is stable to multiple freeze/thaw cycles. Be sure to heat the reagent in a 37°C water bath and mix the reagent to ensure a homogenous solution before use.

Q: Do I need to protect Alamar Blue reagent from light?

A: Yes, Alamar Blue reagent is very slowly converted into a fluorescent product over time, when exposed to light, thus leading to high background values. Store the reagent, protected from light.

Methods Questions

Q: What is the optimal incubation time and temperature of cells with Alamar Blue reagent?

A: Incubate the cells with Alamar Blue reagent for 1–4 hours at 37°C. For more sensitive detection with low cell numbers, increase the incubation time for up to 24 hours.

Q: Can you incubate cells with Alamar Blue reagent overnight?

A: Yes. However, signals from higher cell density samples may have “saturated,” which means the linearity of reagent may have reached a plateau. If this occurs, decrease the incubation time.

Q: What if I don't have an instrument suitable for reading fluorescence?

A: The absorbance of Alamar Blue reagent also changes depending on cell viability and proliferation. Therefore, simply monitor the absorbance of the reagent at 570 nm, while using 600 nm as a reference wavelength.

Q: Is Alamar Blue assay strictly an endpoint assay?

A: No. While Alamar Blue can be used as a terminal readout of a population of cells, the reagent can also be used to continuously monitor cell viability and proliferation in real time. Since Alamar Blue reagent is non-toxic, you can incubate cells with reagent and monitor fluorescence (or absorbance) over time on the same sample.

Troubleshooting Questions

Q: What is the problem for observing high background fluorescence values?

A: The reagent may be breaking down due to exposure to light. Be sure to store Alamar Blue reagent in the dark and do not expose the reagent to direct light for long periods of time.

Q: Why are the fluorescence values so low in intensity?

A: Try increasing the incubation time of cells with Alamar Blue reagent, changing the instrument's “gain” setting, and checking the instrument filter/wavelength settings. Make sure to have positive controls (living cells) in the experimental design for troubleshooting.

Q: Why are the fluorescence values so high that they are beyond the linear range of the instrument?

A: Try decreasing the incubation time or reducing the number of cells used in the experiment.