



全谱绝对定量代谢组学研究

池逸 Ph.D.

ychi@metaboprofile.com

香港浸会大学 贾伟教授课题组



专注代谢组学服务

Q300全谱绝对定量代谢组学

样本类型:

血清/血浆: 20uL


血片: 3个

尿液: 20uL

冻干粪便: 5mg


组织: 10mg

细胞: 10^6



For FTN & FL: 96 well plate (Waters 700 uL round well, Mfr#: 186005837; 7 mm round plug cap-mat, pre-slit, Mfr#: 186006332)

	1	2	3	4	5	6	7	8	9	10	11	12
A 01	S01	S02	S03	S04	S05	S06	S07	S08	S09	S10	S11	S12
B 13	S13	S14	S15	S16	S17	S18	S19	S20	S21	S22	S23	S24
C 25	S25	S26	S27	S28	S29	S30	S31	S32	S33	S34	S35	S36
D 37	S37	S38	S39	S40	S41	S42	S43	S44	S45	S46	S47	S48
E 49	S49	S50	S51	S52	S53	S54	S55	S56	S57	S58	S59	S60
F 61	S61	S62	S63	S64	S65	S66	S67	S68	S69	S70	S71	S72
G 73	S73	S74	S75	S76	S77	S78	S79	S80	S81	S82	S83	S84
H	CAL0	CAL1	CAL2	CAL3	CAL4	CAL5	CAL6	CAL7	QC1	QC2	QC3	QC4
	85	86	87	88	89	90	91	92	93	94	95	96



类别	代谢物	回收率%	CV%	类别	代谢物	回收率%	CV%
氨基酸类	甘氨酸	106	13.3	胆汁酸	胆酸	98	3.3
	丙氨酸	89	2.6		脱氧胆酸	93	3.1
	瓜氨酸	92	2.7		牛磺胆酸	110	9.8
	酪氨酸	103	0.6		牛磺脱氧胆酸	103	4.7
	色氨酸	108	2.3		甘氨酸胆酸	87	2.6
	亮氨酸	102	0.5		甘氨酸脱氧胆酸	92	7.6
脂肪酸类	DHA	93	5.7	还原糖	葡萄糖	84	6.2
	EPA	110	2.8		果糖	80	4.5
	ARA	102	1.2	肉碱类	乙酰肉碱	114	0.8
	油酸	98	4.6		丙酰肉碱	106	0.4
	肉豆蔻酸	116	12.8		戊二酰肉碱	96	1.2
有机酸	柠檬酸	82	1.2	SCFA	丙酸	83	2.2
	琥珀酸	99	0.6		丁酸	87	1.4
	延胡索酸	98	4.5	苯环型化合物	苯乙酸	97	0.6
	2-羟基戊二酸	91	2.2	吡啶类	吡啶乙酸	98	0.9

靶向代谢物组批量定量系统 (TMBQ)

登录

用户名

密码

登录

专利申请号：
CN201910170381.0

物质信息配置表 样本列表 原始数据

任务管理

任务名	任务状态	开始时间	结束时间	操作
project_2020_01_13_13_24_49	成功!	2020-01-13 13:25:58	2020-01-13 14:17:59	🔍 🗑️ 🔄
project_2020_01_08_11_52_18	成功!	2020-01-08 11:53:43	2020-01-08 12:48:08	🔍 🗑️ 🔄

显示第 1 到第 2 条记录, 总共 2 条记录

sample	batch	Pyro	Hydro	Pico	Ortho	Hyc	Hippo	N	Acetyl	L	Malic	εN	Acetyl	N	Acetyl	Indole	3
20200211_q300_std_0	1	0.254171	0.606629	0.194784	0.22616	0.432359	0	0.272379	0.00688	0.106196	0.112958	0.19146					
20200211_q300_std_1	1	0.371957	0.424853	0.326503	0.344752	0.522735	0	1.646386	0.056131	0.17634	0.148425	0.319297					
20200211_q300_std_2	1	0.442377	0.376111	0.431891	0.479824	0.577133	0.048738	1.29509	0.294604	0.19796	0.097184	0.38475					
20200211_q300_std_3	1	1.642246	1.306746	1.412678	1.506927	1.471589	1.108071	1.506797	1.182362	1.329138	1.598956	1.394449					
20200211_q300_std_4	1	2.940201	3.360522	2.932135	2.896044	2.882367	2.866398	2.258871	2.750231	2.723217	2.817188	3.014732					
20200211_q300_std_5	1	5.717008	5.811097	6.08098	6.096581	5.701514	6.302752	4.294303	6.395314	6.149966	6.982818	6.114051					
20200211_q300_std_6	1	12.32225	11.5123	12.42311	12.10964	11.90451	12.75175	6.948525	12.71619	12.63291	12.1511	12.37478					
20200211_q300_std_7	1	25.21603	25.50799	25.10417	25.24602	25.41374	27.54333	11.01095	26.16887	26.49371	24.99763	25.11273					
20200211_q300_std_8	1	46.84758	45.23579	47.12114	44.11559	45.38827	48.50795	50.13242	49.33566	49.20301	54.62551	44.43227					
095_126_qc_01	1	0.283921	3.018995	0.610209	4.550714	2.534365	270.8913	30.25411	4.475557	43.94139	6.681148	0.537864					
095_126_qc_02	1	0.267985	2.498426	0.563047	4.094282	2.643134	275.6561	22.19272	4.017863	41.7003	6.276529	0.59341					
095_126_qc_03	1	0.269749	2.724525	0.619172	3.981024	2.646183	279.5116	0.127958	4.078977	43.88621	6.221438	0.606274					
095_126_qc_04	1	0.339415	2.761152	0.645408	4.130333	2.713694	283.961	19.76919	4.374255	47.17243	6.749566	0.70109					
095_126_qc_05	1	0.324016	3.158008	0.62772	3.791707	2.813336	301.4185	24.26018	4.498699	50.06503	5.385426	0.722045					
095_126_qc_06_1	1	0.318187	3.125866	0.603568	3.345127	2.562364	280.0124	12.6236	4.321499	47.41127	4.562225	0.76444					
095_126_qc_06_2	1	0.298317	3.176583	0.618806	2.952006	2.727718	272.2842	0.127958	4.276905	47.28627	3.981315	0.66249					
095_126_qc_07	1	0.312141	2.93003	0.594425	2.717806	2.59096	272.5818	15.38081	3.989414	44.26271	4.384391	0.70132					
202087	1	0.454913	0.742335	1.008318	9.968044	2.715577	417.1332	0.127958	1.551622	47.7034	3.179827	0.677291					
454619	1	0.281521	1.010324	0.851437	2.253026	2.003447	273.1357	5.837075	3.513766	35.96467	4.235056	0.593492					
579478	1	0.278348	1.841455	0.724646	5.061829	1.582465	299.0768	9.014818	2.468872	27.83098	3.839679	0.634436					
666145	1	0.27654	1.266798	0.341723	1.962915	1.317555	283.6014	10.7077	1.75665	26.09057	2.830612	0.818154					
674924	1	0.378524	0.089374	0.428963	1.628079	0.92232	85.60133	2.946531	0.49028	8.324565	2.362477	0.653083					
803410	1	0.282337	8.328133	0.484397	0.713195	1.072852	67.7568	1.603667	10.00563	26.77904	1.524139	0.173834					
805051	1	0.30236	0.422016	0.780717	5.247738	2.576321	325.2342	0.127958	4.392995	48.56887	3.453319	0.405382					
809405	1	0.27598	0.454917	0.467932	1.199075	2.477845	252.9896	0.127958	6.521198	49.95894	3.895528	0.558597					
958587	1	0.28601	0.972789	0.676851	2.29626	3.425966	239.7322	1.833958	12.45146	94.20323	4.816237	0.743212					
983192	1	0.328072	1.265286	0.712124	5.919883	5.722916	500.7877	13.50295	0.474762	74.89486	7.503691	3.089653					
989283	1	0.253256	0.016109	0.355573	1.52284	1.202401	160.255	0.127958	0.899524	5.185305	2.036943	0.369949					
989439	1	0.269623	0.279681	0.298223	0.975956	1.064245	30.87071	0.127958	0.711578	7.383481	1.058485	0.353336					
989484	1	0.272511	0.668104	0.370871	1.246964	1.394101	214.7472	0.127958	1.868306	16.63349	2.332363	0.362609					
991303	1	0.261712	0.128887	0.428224	2.020712	1.537912	287.5538	19.12065	3.035148	31.5919	4.737505	0.569615					
991378	1	0.292276	0.54071	0.362719	2.303814	1.713669	283.1211	3.174366	0.628647	18.19773	3.099177	0.608094					
991550	1	0.281038	0.836987	0.345088	1.284438	1.221288	224.7506	8.440007	5.06714	11.75053	3.253181	0.599553					







图1 多变量质控图




✓ 3000+例血浆队列样本




✓ 1000+例粪便队列样本

类别	血清/血浆	粪便	尿液	组织	干血斑
氨基酸	44	42	46	45	35
有机酸	27	32	32	32	29
脂肪酸	45	34	24	46	35
胆汁酸	17	32	18	20	5
碳水化合物	20	17	17	16	17
苯环类	17	13	16	10	8
肉碱	21	13	13	18	19
苯丙酸类	4	8	6	/	/
短链脂肪酸	8	8	5	5	6
吡啶类	3	2	4	/	/
其他	3	4	19	15	11
总代谢物	209	205	200	207	165


检测



定量



分析



靶向代谢物组批量定量系统 (TMBQ)

登录

用户名

密码

登录

云平台iMAP

- 多维统计
- 单维统计
- 通路分析
- 诊断分析
- 关联分析

Waters
THE SCIENCE OF WHAT'S POSSIBLE.™

AB SCIEX

ThermoFisher
SCIENTIFIC

Agilent Technologies

岛津
SHIMADZU







ВУНКЕВ

MP 麦特宏谱
Metabo-Profile




目录


- 目录 1
- 一、项目简介 7
- 二、质量控制 8
 - 1 多变量质控图 8
 - 2 QC样本间相关性 9
 - 3 样本代谢轮廓变异 10
 - 4 代谢物鉴定及注释情况 11
- 三、分析结果 19



Gut Microbiota




PubMed文章数量



2017年首次建立基于GC-TOFMS的自动化全定量肠道菌群代谢物分析方法

High Throughput and Quantitative Measurement of Microbial Metabolome by Gas Chromatography/Mass Spectrometry Using Automated Alkyl Chloroformate Derivatization

Linjing Zhao,^{†,‡,#} Yan Ni,^{†,§,#} Mingming Su,^{†,§,#} Hongsen Li,[§] Fangcong Dong,[§] Wenlian Chen,[§] Runmin Wei,[§] Lulu Zhang,[§] Seu Ping Guiraud,^{||} Francois-Pierre Martin,^{||} Cynthia Rajani,[§] Guoxiang Xie,[§] and Wei Jia^{*,†,§}



- 短链脂肪酸 (12)
- 胆汁酸 (40+)
- 苯类 (20+)
- 吡啶类 (10+)
- 苯丙酸类 (7)
- 氨基酸 (40+)
- 中长链脂肪酸 (40+)
- 碳水化合物 (10+)
- 有机酸 (30+)
- 其他 (5)

Gut microbiota–bile acid–interleukin-22 axis orchestrates polycystic ovary syndrome

Xinyu Qi^{1,11}, Chuyu Yun^{1,2,11}, Lulu Sun^{1,2}, Jialin Xia^{1,2}, Qing Wu^{1,2}, Ying Wang^{1,3,4}, Lina Wang^{1,3,4}, Yangming Zhang^{1,2}, Xianyi Liang^{1,2}, Liying Wang^{1,3,4}, Frank J. Gonzalez⁵, Andrew D. Patterson⁶, Huiying Liu^{1,2}, Liangshan Mu^{1,3}, Zehong Zhou^{1,3}, Yue Zhao^{1,3,7}, Rong Li^{1,3,4,7}, Ping Liu^{1,3,4,7}, Chao Zhong⁸, Yanli Pang^{1,3,7*}, Changtao Jiang^{1,2*} and Jie Qiao^{1,3,4,7,9,10*}

Ma *et al.*, *Science* **360**, eaan5931 (2018)

Gut microbiome–mediated bile acid metabolism regulates liver cancer via NKT cells

Chi Ma, Miaojun Han, Bernd Heinrich, Qiong Fu, Qianfei Zhang, Milan Sandhu, David Agdashian, Masaki Terabe, Jay A. Berzofsky, Valerie Fako, Thomas Ritz, Thomas Longerich, Casey M. Theriot, John A. McCulloch, Soumen Roy, Wuxing Yuan, Vishal Thovarai, Shurjo K. Sen, Mathuro Ruchirawat, Firouzeh Korangy, Xin Wei Wang, Giorgio Trinchieri, Tim F. Greten*

初级胆酸可激活CXL蛋白，进而激活NKT细胞，抑制肿瘤，而次级胆汁酸有相反的作用。

GDCA通过激活肠道3型固有淋巴细胞（ILC3）的GATA3通路刺激IL-22分泌，进一步促进白色脂肪棕色化以及抑制卵巢局部炎症，进而改善PCOS样表型。

Nature, 2019, doi:10.1038/s41586-019-1785-z

Article

Bile acid metabolites control T_H17 and T_{reg} cell differentiation

3-oxo-LCA可抑制Th17细胞的分化，而isoalloLCA可增强Treg细胞的分化。

Alzheimer's
&
Dementia

Alzheimer's & Dementia 15 (2019) 76-92

Alzheimer's & Dementia ■ (2018) 1-13

Altered bile acid profile in mild cognitive impairment and Alzheimer's disease: Relationship to neuroimaging and CSF biomarkers


Kwangsik Nho^{a,1}, Alexandra Kueider-Paisley^{b,a}, Siamak MahmoudianDehkordi^{b,a}, Matthias Arnold^{b,c}, Shannon L. Risacher^a, Gregory Louie^b, Colette Blach^d, Rebecca Baillie^c, Xianlin Han^f, Gabi Kastenmüller^{c,g}, Wei Jia^h, Guoxiang Xie^h, Shahzad Ahmadⁱ, Thomas Hankemeier^j, Cornelia M. van Duijnⁱ, John Q. Trojanowski^k, Leslie M. Shaw^k, Michael W. Weiner^l, P. Murali Doraiswamy^{b,m,n}, Andrew J. Saykin^{a,**}, Rima Kaddurah-Daouk^{b,m,n,*}, for the Alzheimer's Disease Neuroimaging Initiative and the Alzheimer Disease Metabolomics Consortium²

胆汁酸可能通过肠-肝-脑轴影响AD的进展。

Altered bile acid profile associates with cognitive impairment in Alzheimer's disease—An emerging role for gut microbiome

Siamak MahmoudianDehkordi^{a,1}, Matthias Arnold^{a,b,1}, Kwangsik Nho^{c,1}, Shahzad Ahmad^d, Wei Jia^{e,f}, Guoxiang Xie^e, Gregory Louie^a, Alexandra Kueider-Paisley^a, M. Arthur Moseley^g, J. Will Thompson^g, Lisa St John Williams^g, Jessica D. Tenenbaum^h, Colette Blachⁱ, Rebecca Baillie^j, Xianlin Han^k, Sudeepa Bhattacharyya^l, Jon B. Toledo^m, Simon Schaffnerⁿ, Sebastian Kleinⁿ, Therese Koalⁿ, Shannon L. Risacher^c, Mitchel Allan Kling^o, Alison Motsinger-Reif^p, Daniel M. Rotroff^p, John Jack^p, Thomas Hankemeier^q, David A. Bennett^r, Philip L. De Jager^s, John Q. Trojanowski^t, Leslie M. Shaw^t, Michael W. Weiner^u, P. Murali Doraiswamy^{a,v,w}, Cornelia M. van Duijn^d, Andrew J. Saykin^{c,**}, Gabi Kastenmüller^{b,x,*}, Rima Kaddurah-Daouk^{a,v,w,***}, for the Alzheimer's Disease Neuroimaging Initiative and the Alzheimer Disease Metabolomics Consortium


▶▶ BAP™ kit (Bile Acid Profiling)介绍



m/z 375	m/z 482	m/z 389
m/z 391	m/z 498	m/z 405
m/z 407	m/z 514	m/z 403
m/z 432	m/z 508	m/z 377
m/z 448	m/z 458	m/z 393
m/z 464	m/z 373	m/z 401
U01 dehydroLCA	U14 HDCA	G05 GHCA
U02 LCA	U15 UDCA	U28 βCA
U03 isoDCA	U16 CA	G06 GUDCA
U04 isoLCA	T02 TDCA	U29 ωMCA
U05 alloLCA	U17 ACA	T04 TCA
G01 GLCA	U18 βHDCA	G07 GHCA
U06 DCA	U19 3-DHCA	T05 THDCA
U07 CDCA	T03 TCDC	T06 TUDCA
T01 TLCA	U20 βUDCA	U30 7,12-diketoLCA
U08 apoCA	U21 muroCA	U31 DHCA
U09 NorDCA	U22 HCA	U32 UCA
U10 12-ketoLCA	U23 12-DHCA	T07 THCA
U11 7-ketoLCA	U24 βMCA	T08 TβMCA
U12 6,7-diketoLCA	U25 NorCA	T09 TαMCA
G02 GDCA	U26 7-ketoDCA	G08 GDHCA
U13 6-ketoLCA	G04 GCA	T10 TωMCA
G03 GCDC	U27 αMCA	T11 TDHCA



- Anal Chem, 2016, 88(14): 7041-7048
- Anal Bioanal Chem, 2017, 409(23):5533-45
- Anal Bioanal Chem, 2018, 410(21):5287-300
- Nat Commun. 2019. 10:4971
- Alzheimer's & Dementia, 2019, 15: 76-92
- Nat Rev Gastroenterol Hepatol, 2018, 15(2):111-28
- Alzheimer's & Dementia, 2018, doi:10.1016/j.jalz.2018.08.012
- EBioMedicine, 2018, doi:10.1016/j.ebiom.2018.10.030
- BMC Biology , 2017,15:120 -34
- International Journal of Cancer, 2016, 139(8):1764-75





非靶向及定量代谢组学

Q300全定量代谢组
宏代谢组
胆汁酸
氨基酸
有机酸
短链脂肪酸
中长链脂肪酸
氧化三甲胺(TMAO)

核苷酸
辅酶
糖代谢
激素
脂质
维生素
还原糖
神经递质



同位素代谢流分析

代谢流分析
同位素示踪分析



微生物组

16S rDNA测序
宏基因组



定量方法检测试剂盒

Q300全定量代谢芯片
Q200宏代谢
胆汁酸
脂肪酸



麦特绘谱部分合作客户/文章

Cell Metabolism, 2020, 31, 1-15 (IF=22.41)

Immunity, 2019, 50, 1-16 (IF=19.73)

Gut, 2019, doi:10.1136/gutjnl-2018-317609 (IF=17.02)

Hepatology, 2020, doi:10.1002/hep.31280 (IF=14.97)

Signal Transduction and Targeted Therapy, 2020, 5: 177-90 (IF=13.49)

Microbiome, 2020, 8: 74 (IF=11.61)

Microbiome, 2020, 8: 133 (IF=11.61)

The EMBO Journal, 2020, 39: e103304 (IF=9.89)

Arthritis Rheumatol, 2020, doi:10.1002/art.41419 (IF=9.59)

Chemical Science, 2020, 11: 6111 (IF=9.35)


Theranostics, 2019; 9(3): 900-919 (IF=8.71)

Theranostics, 2019; 9(10): 2999-3013 (IF=8.71)

Clin Gastroenterol Hepatol, doi:10.1016/j.cgh.2020.06.067 (IF=8.55)

Carbohydrate Polymers, 2020, 248, 116780-92 (IF=7.18)

Oncogene, 2020, 39, 2358-76 (IF=6.64)



优质服务

不断超越客户期望
建立长期合作关系

专业化人才队伍
不断创新,追求新高度

精益求精



科学管理

全方位科学管理
为客户提供专业产品



Thank you